Electronic Lab Notebook (ELN) Entry

Technician Name: Emily Saunders

Date: 2025-03-15

Batch ID: WHT-HRS-4522

Sample Type: Hard Red Spring Wheat

Analysis Performed: Moisture & Protein Content

Instruments Used: Analytical Balance (ML204), Oven (UN55), NIR Analyzer (FOSS 1241)

Start Time: 09:15 AM

End Time: 11:40 AM

Observations: Sample was finely ground before analysis. Moisture oven temperature held

stable at 105°C throughout.

Initial sample weight: 5.032g | Final weight after drying: 4.410g

Moisture calculated: 12.35% (within spec).

Protein reading via NIR: 13.6% (well above minimum threshold).

Notes: Results consistent with vendor report. No deviations observed. All readings

logged in ELN and verified by technician.

Reviewed By: Mike Skelley, QA Analyst

Review Date: 2025-03-15